



Operation Manual



H & R Healthcare
1750 Oak Street
Lakewood, NJ 08701
Phone: 732-367-5533
Fax: 732-367-6231
www.handrhealthcare.com



Congratulations and thank you for purchasing this Signa Relief Alternating Pressure System with Low Air Loss. PLEASE READ THIS OPERATION MANUAL CAREFULLY BEFORE SETTING UP AND USING THE DEVICE. Pay special attention to the warnings and other safety information. Use of genuine manufacturer components is essential for optimal performance. If you do not fully understand all the instructions, safety precautions, and warnings, do not use this device. If you have any questions or concerns, please contact H & R Healthcare at 732-367-5533.

Table of Contents

Section	Description	Page
1.0	Indications	4
2.0	Safety Precautions.....	4
3.0	Cautions	5
4.0	System Package.....	5
5.0	Features.....	6
6.0	Mattress and Power Unit Installations.....	7
7.0	Program Settings.....	8
8.0	Patient Transfer and Transport.....	8
9.0	Emergency CPR Deflation.....	9
10.0	Cleaning Instructions.....	9
11.0	Routine Maintenance.....	9
12.0	Troubleshooting.....	10
13.0	Returns for Service	11
14.0	Warranty.....	11
15.0	Product Specifications.....	12
16.0	EMC Declaration	13

1.0 Indications



Operation Instructions



Caution



Declaration of Conformity to Medical Device Directive.

IP21

Do not immerse power unit in liquid or spray liquids directly on power unit.



Type BF Applied part



Protective Earth (PE)



NRTL_SGS Q Mark Listing logo



WEEE Logo

Indications

This alternating pressure system with low air loss is designed to treat and prevent wounds by facilitating blood circulation and decreasing pressure of each tissue's contact area. Always consult a physician or health professional before using this mattress system.

Contraindications

This device is not indicated for treatment of certain patient conditions such as fracture of unstable vertebrae and illness of unstable vertebrae. Always consult a physician or health professional before using this device.

The use of this system does not replace the regular repositioning, monitoring, and nursing of the patient.

2.0 Safety Precautions

1. To ensure proper operation, please inspect and verify that all parts are set up properly and are anchored securely. Do not place anything on top of the power unit. Make sure power cord is underneath bed frame and does not pose a hazard.
2. It is recommended to limit bed linens to a single layer in order to allow moisture to escape efficiently through the coverlet. Breathable incontinent pads are recommended for use with this mattress system.
3. Avoid using this device near open flames, lighters, or cigarettes. Flammability hazard exists. This device draws air from the surrounding environment. Thus, cigarette smoking may damage internal components.
4. This system should be disinfected thoroughly between patients to avoid cross contamination.
5. Verify that patient weight does not exceed the weight capacity of the bed frame, bed rails, or this mattress system.

3.0 Cautions

1. *Use this mattress with proper side rails to ensure that the gap between the side rail and the top of the mattress is not large enough to pose risk of head or neck entrapment. Failure to do so could result in serious patient injury.*
2. DO NOT disassemble the power unit if you are not a qualified technician. Please contact H & R Healthcare for all service.
3. This product is NOT AP/APG protected.
4. It is recommended that the patient be repositioned periodically while using this mattress.
5. Incorrect battery replacement may result in risk of explosion. Replace only with same or equivalent battery type recommended by the manufacturer.
6. If the equipment needs maintenance, contact H & R Healthcare as soon as possible. For equipment that is no longer functional, make sure to follow national and/or state requirements for disposal of the unit.
7. Ensure that any disposal of the unit, accessories, waste products, residues etc. complies with national and local law.

8. The power plug serves to disconnect the device and should not be positioned in a way that makes it difficult to disconnect.
9. No modification of this equipment is allowed.
10. To avoid the risk of electric shock, this equipment must be connected to a supply mains with protective earth.
11. To shut down the unit, press the power button.

CAUTION: ENSURE THAT THERE ARE NO PROTRUDING OBJECTS, SHARP POINTS, OR BED SPRINGS UNDER THE MATTRESS AS THESE COULD PUNCTURE THE AIR CELLS.

4.0 System Package

Power Unit Package

- Power Unit x 1
- Power Cord Set x 1
- Operational Manual x 1


Mattress Package

- Mattress replacement unit with coverlet x 1

Transport & Package

- Temperature: -25°C to +70°C
- Relative Humidity: 0% to 90%

5.0 Features



Control Panel Features

Power/Mute

The green power LED will light when the main power rocker switch on the right side panel is turned on.

Weight (lbs)

Simply press weight button to adjust patient weight from 100 lbs to 325 lbs according to patient weight. The scale is only an approximation. Please adjust the weight setting if the mattress is too soft or firm to suit each patient's needs.

Caregivers should always perform a hand check by placing their hands underneath patient's pelvis area to check if there is sufficient air

support to ensure the patient is not bottoming out.



Alert LED

This alert notifies medical staff the mattress has insufficient pressure. Audible alarm will sound for 15 seconds at each cycle, and the indicator will stay on until the problem has been solved.




Pressure monitoring LED

Power unit features an integrated pressure sensor which monitors the mattress' internal pressure to ensure maximum pressure relief. LED indicates that the system is monitoring the mattress' internal pressure while the compressor is inflating.

Alternating/Static

Press Mode to change between alternating and static mode. Alternating mode will have a cycle time of 10 minutes.

Right Side Panel Features



Rocker Switch (1)

Power unit's main power switch.


Couplers (2)

Quick release female couplers are used to secure mattress air hoses to power unit.

Power Receptacle (3)

Insert power cord firmly into receptacle.

Rear Panel Features



Convenient Handle (1)

The handle provides additional gripping surface for user to carry the power unit.

Hanging Hooks (2)

Hanging hooks are designed to hang the power unit on almost any foot board.

Fuse (3)

Fuse Holder

Air Filter and Cap (4)

H & R Healthcare recommends that the filter be kept clean to ensure optimal performance of the power unit.

6.0 Mattress and Power Unit Installations

1. Remove existing mattress from bed frame.
2. Place H & R Healthcare mattress replacement on the bed frame with the logo at foot end. Secure the mattress at each side using anchor straps. Ensure that the anchor straps do not interfere with any moving parts of the bed frame before proceeding to the next step.
3. Secure the power unit onto the foot board using hanging hooks.
4. Ensure that the mattress; CPR latch is firmly inserted and the CPR tag is in place.

5. Firmly connect the air hose couplers to the couplers on the power unit.
6. Plug the power unit into an electrical outlet, and turn on the main power rocker switch on the right side of the power unit.
7. Press the power/mute button to inflate the mattress system. Please wait for approximately 30 seconds until the power unit is ready before making any program adjustments.
8. Wait approximately 30 minutes for the mattress to inflate fully.

7.0 Program Settings

1. Place the patient in the center of the mattress. Adjust the mattress' internal pressure according to the patient weight by using the weight button on the control panel of the power unit. If the mattress is too soft or too firm, increase or decrease the mattress' internal pressure one increment at a time and wait for the system to stabilize before making another change. Continue this process until comfort is achieved.
2. Caregivers should always perform a hand check by placing their hands underneath patient's pelvis area to check if there is sufficient air support to ensure the patient is not bottoming out.

3. **IMPORTANT: TUCKING THE SHEET IN TOO TIGHTLY REDUCES THE EFFECTIVENESS OF THE SYSTEM.**

8.0 Patient Transfer and Transport

Transfer

It is recommended to have the mattress system at maximum weight setting during transfer. Make sure the bed is secured before proceeding.

Transport

In the event of a patient transport, two options are available.

1. Detach the mattress' air hose couplers from the power unit's quick release couplers and connect the two air hose couplers together to retain air in the mattress. The mattress will stay inflated for approximately 2 hours, depending on the patient's weight.
2. Unplug the power unit's power cord from the wall outlet and the mattress should stay inflated for approximately 40 - 50 minutes. This mattress system also has a safety foam base underneath the air cells to support the patient for a short period of time in case of deflation. To resume normal operation, please follow 7.0 Program Settings.

9.0 Emergency CPR Deflation

In case of emergency, pull hard to remove the mattress' CPR latch for emergency deflation, and turn off the power unit by pressing the power/mute button on the control panel at the same time. The mattress will deflate from the weight of the patient. To resume normal operation, simply reinsert the CPR latch securely, press the power/mute button again and reset the patient weight.

10.0 Cleaning Instruction

Air mattress and power unit must be cleaned thoroughly between patients to avoid cross contamination. The following is a suggested guideline. Be sure to follow local infection control policies.

Regular cleaning can be performed at bedside with disinfectant followed by drying with a clean dry cloth. Use only mild detergents to clean the coverlet and the mattress. Any appropriate NON-PHENOLIC cleaning agent may be used for heavy soiling from urine, blood, or

other bodily fluids. Please ensure that air mattress and coverlet are completely dry before letting the patient lie on the surface again.

Do not use electric or tumble dryers. Do not iron.

WARNING! Always unplug the power unit before cleaning. Routine cleaning of power unit can be done by wiping down with damp cloth using disinfectant or mild detergent. Never spray liquids directly on the unit itself.

11.0 Routine Maintenance

Remove air filter from the rear panel of the power unit by opening up the filter cap. Inspect the filter for dirt or dust, and clean it with mild soap and water. Reinsert the dried air filter after cleaning and ensure that the cap is secure. If replacement is needed, contact H & R Healthcare.

Only disinfected and dry systems should be stored. Disconnect the air hoses from the power unit. Roll up the mattress starting from the head end and working down toward the foot end. Use the straps to secure it.

12.0 Troubleshooting

Problem	Inspection Procedure	Possible Solutions
1. Power unit is not working.	Check if power cord is firmly plugged into both the control unit and the electrical outlet.	Secure power cord into control unit and/or electrical outlet.
	Check if the power switch is in the ON position.	Turn power switch to ON position.
	Check if power surge has shut down the power unit.	A power surge may overload the circuitry temporarily. Turn the unit off, and check the fuse for damage. Turn the unit on again with normal procedure.
	Make sure there is no power failure.	
	Power unit does not respond to possible solutions.	Please contact H & R Healthcare for assistance.
2. Low Pressure alert does not turn off after 30 seconds after the system is turned on.	Possible control failure.	Please contact H & R Healthcare for assistance.
3. Alarm LED is lit during operation.	Check if there is leakage in air tubes or air cells	Please contact H & R Healthcare for assistance.
4. Power unit is working, but mattress replacement is not inflating and/or Bottoming out is occurring and/or Patient leaves a deep indentation at the contact area which does not return back to its original shape.	Verify that CPR latch is properly inserted.	Insert CPR latch securely.
	Check if mattress' air hose couplers are properly connected to power unit's quick release couplers.	Secure air hose couplers firmly into place.
	Verify that patient weight setting is correct.	Increase or decrease weight setting until appropriate pressure is reached.
	Inspect air filter for dust.	Clean or replace air filter.

	Lift mattress coverlet up to check if air cells are connected correctly.	Make sure all air cells are properly linked to air supply.
	Lift mattress coverlet up to check if air tubes are kinked or obstructed.	Check and adjust air tubes' positions.
	Check if air cells are cut or cracked.	Please contact H & R Healthcare for assistance.
5. Patient's wounds are not responding to pressure relief (reddening of skin).		Contact your physician and/or nursing service immediately.

13.0 Returns for Service

Service and repair must be performed by an H & R Healthcare authorized technician or representative. Please contact H & R Healthcare at 732-367-5533.

14.0 Warranty

H & R Healthcare warrants this product against manufacturer defects for one year from the date of placement.

Please inspect this product upon purchase and notify H & R Healthcare if there is any damage or missing parts.

Warranty coverage will not be extended to any product on which the production lot number has been removed or defaced, on which repair has been attempted by any person or agency not authorized by H & R Healthcare or, if in the sole

opinion of H & R Healthcare, the system shows evidence of tampering, abnormal or unreasonable abuse, negligence, accident or operation without regard for the restrictions specified in the instructions which accompany the system. This warranty does not cover routine maintenance such as cleaning, adjustment, lubrications, and updating of equipment or parts. If, in the sole opinion of H & R Healthcare, the damage is not covered under this warranty, the purchaser will be responsible for labor and parts costs.

The warranty stated above is the only warranty made and is in lieu of all other warranties whether expressed or implied. H & R Healthcare, nor any of its agents, shall be held liable for consequential or incidental damages of any kind.

15.0 Product Specifications

Air Pump (Power Unit) Model
SR366

Mattress Dimension
80"/84" x 37" x 7"

Number of Air Cells
20 Air Cells

Power Unit Dimension
11.5" x 8.5" x 4.5"

Integrated Foam
2"

Material
Power Unit - Plastic Case rated UL 94V-0
Mattress Cell - 100 % nylon with TPU
lamination
Coverlet - 100% nylon quilted
backing

Power Input
100-120Vac 60Hz 0.2A

Power Consumption
Normal Operation: Max 15W

Power Control Unit Weight
5.5 lbs

Fuse Rating
T1A 250V

Electrical Classification
Class I Type BF. Applied part is
mattress. Continuous operation.

IP21, do not immerse power unit in
any liquid or spray any liquids
directly on the power unit.
This system is not AP/APG
Protected.

Safety Test & Certified Standard
NRTL_SGS Q Mark Listing to
ANSI/AAMI ES60601-1:2005+A1:2012
CAN/CSA C22.2 No.60601-1:14
EN 60601-1:2006/A1:2013_v3.1
IEC 60601-1:2005/A1:2012_v3.1

EMC Test Standard
IEC 60601-1-2:2014/EN60601-1-2:2015
Emission Classified as Class B to CISPR 11

Operating Conditions
Temperature range:5°C to 40°C
Relative humidity range: 15% to 93%
Atmosphere range: 700hPa - 1060hPa

Maximum Weight Capacity
325 lbs

**H & R Healthcare reserves the
right to change prices or design
parameters and specifications
without incurring any liability in
connection with such
modifications.**

16.0 EMC Declaration

Declaration of Conformity


For EN 60601-1-2

Company Name:	CARILEX MEDICAL, INC
Company Address:	No. 77, KEJI 1ST RD. GUSHAN DIST., TAOYUAN CITY,333 TAIWAN (R.O.C.)
Trade Name:	Carilex
Report Number:	ETC 16-08-RBO-002-01
Power Supply:	AC 100 - 120V 60Hz 0.2A 20W max.

Recommended separation distances between portable and mobile RF communications equipment and the ME equipment			
The Air Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Air Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Air Pump as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

Declaration – electromagnetic emissions		
The Air Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Air Pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Air Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

**Declaration – electromagnetic emissions and immunity –
for EQUIPMENT and SYSTEMS that are use in the professional healthcare facility
environment or in the home healthcare environment**

The Air Pump declaration – electromagnetic immunity					
The Air Pump system is intended for use in the electromagnetic environment specified below.					
The customer or the user of the Air Pump system should assure that it is used in such an environment.					
Immunity test	IEC 60601 test level		Compliance level		Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms ; 6 Vrms 150 kHz to 80 MHz		3 Vrms ; 6 Vrms 150 kHz to 80 MHz		Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%		3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%		
Proximity fields from RF wireless Communications equipment IEC 61000-4-3	27 V/m	385 MHz	27 V/m	385 MHz	Interference may occur in the vicinity of equipment marked with the following symbol. 
	28 V/m	450 MHz	28 V/m	450 MHz	
	9 V/m	710 MHz	9 V/m	710 MHz	
		745 MHz		745 MHz	
		780 MHz		780 MHz	
	28 V/m	810 MHz	28 V/m	810 MHz	
		870 MHz		870 MHz	
		930 MHz		930 MHz	
	28 V/m	1720 MHz	28 V/m	1720 MHz	
		1845 MHz		1845 MHz	
		1970 MHz		1970 MHz	
	28 V/m	2450 MHz	28 V/m	2450 MHz	
9 V/m	5240 MHz	9 V/m	5240 MHz		
	5500 MHz		5500 MHz		
	5785 MHz		5785 MHz		

Declaration – electromagnetic immunity			
The Air Pump system is intended for use in the electromagnetic environment specified below.			
The customer or the user of the Air Pump system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV , ±4 kV , ±8 kV , ±15 kV air	±8 kV contact ±2 kV , ±4 kV , ±8 kV , ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV ±1 kV differential mode ±2 kV common mode	±0.5 kV ±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T ; 0 , 5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°	0 % U_T ; 0 , 5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.