


Quick Reference Guide


– Convertible Chairs

Description

Our convertible chair is an unique device that easily converts from a stretcher to an upright mobile chair in seconds. Its Tilt-in-Space capability allows for increased blood circulation helping to prevent skin breakdown and other complications of immobility. It provides a no-lift solution for moving, repositioning, and transferring patients in the safest possible manner.



Transfer a patient: TIP: Always lock wheels on chair before attempting a transfer.



Tilting (down or up):

<p>Hold the headrest handle with one hand and squeeze the tilting lever and push handle together with other. Push both slowly down or lift slowly upward. (keeps patient in the desired seated angle)</p>	<p>Release the tilting lever when the tilting reaches the desire angle.</p>	<p>TIP: It is recommended to maintain holding handle until after you release any levers. Always check that chair has locked into the position that has been set.</p>	

Quick Reference Guide

– Convertible Chairs

Converting (Chair or fully reclined):

<p>Hold the headrest handle with one hand and squeeze the reclining lever and push handle together with other. Push both slowly down or lift slowly upward to desired position.</p>	<p>TIP: It is recommended to put chair in full upright position (no tilt) before reclining unless the trendelenburg position is desired.</p>

Caster Control Pedal Function:

<p>All Wheel Lock – locks all wheels in place. Pedal tilted away from chair. (check that the chair is locked after applying brake)</p>	<p>Free Wheel Mode – allows wheels to move independently. Pedal is flat.</p>	<p>Steering Mode – locks two front wheels in position to facilitate steering. Pedal is tilted towards chair.</p>

Headrest Support:

<p>Headrest depth and angle can be adjusted by turning levers. This allows for headrest to telescope in and out.</p>	<p>Hold the middle part of the headrest support with one hand and turn one side at a time slowly inward with the other hand.</p>

Adjusting Levers:

<p>Function 1 Tighten/Loosen</p> <p>Function 2 Adjust lever position</p>
<p>Headrest height, headrest depth/angle, footrest height are three levers in the same form with two functions: Function 1 - Tighten and loosen (clockwise to tighten) Function 2 - Change lever position in a tightened state.</p>

99518-Quick reference guide-Convertible Chairs-Rev02 (2015-05-18)



IMPORTANT! This quick reference guide does not replace the user manual which can be downloaded at www.humancaregroup.com

H&R HEALTHCARE

Support Surfaces and Negative Pressure Wound Therapy

1750 Oak Street • Lakewood, NJ 08701 • (800) 801-5533 • www.handrhealthcare.com