



Broda chairs provide exceptional comfort and ergonomic design for individuals requiring long term care. Broda chairs allow for safe and comfortable positioning without the use of restraints.

SASHAY PEDAL CHAIR

- Designed to encourage self mobility
- Up to 18° of infinitely adjustable seat tilt
- Front pivot seat tilt helps maintain the user's foot on floor position for safe self-propelling while tilted

LATITUDE PEDAL ROCKER

- Gentle rocking provides calming, soothing effect for agitated patients
- Up to 18° degrees of infinitely adjustable seat tilt



Latitude Pedal Rocker

SYNTHESIS TILT RECLINER

- Broda's full-featured chair providing exceptional comfort and repositioning
- Up to 40° of infinitely adjustable seat tilt
- Up to 90° of infinitely adjustable back recline
- Trendelenburg positioning



Synthesis Tilt Recliner



Comfort Tension Seating®

- Redistributes pressure to help prevent skin breakdown
- Creates air flow for heat and moisture exchange
- Contours to the body for long term sitting comfort
- Memory retention allows straps to return to their original shape within seconds for long term and multi resident use

Additional models and accessories available. Swing-away leg rests included.