



H E A L T H C A R E



Compliance. Safety. Function.



At its lowest setting, the FloorBed® is just 2.8" tall, minimizing the risk of serious injury if users fall out of bed. The FloorBed® also rises to full nursing height for care duties and is compliant with the latest safety standards.

- **Improves Quality of Life**
No need for bed rails or crash mats
- **Helps Balance Workload**
Lower staff monitoring time
- **Reduces Costs**
Nearly eliminates the risk of bed fall injuries



Mattress platform height
of just 2.8 inches



Backrest and knee brake
functions



Full nursing height,
No trailing wires

FloorBed® 1 Technical Specifications

Dimensions	36.4" x 90.4"	Lifting Range	2.8" - 25.6" (excl. mattress)
Mattress Platform Size	Standard configuration: 36.0" x 80.0" Extended configuration*: 36.0" x 84.0"	Profiling functions	Backrest with Auto-regression, Knee brake, Trendelenburg, Reverse Trendelenburg
Safe Working Load	408 lbs	Max Patient Weight	330 lbs

*Mattress platform extension infill piece MUST be used if the mattress platform is extended.

FloorBed® 1 Plus Technical Specifications

Dimensions	42.3" x 90.4"	Lifting Range	3.1" - 25.6" (excl. mattress)
Mattress Platform Size	Standard configuration: 42.0" x 80.0" Extended configuration*: 42.0" x 84.0"	Profiling functions	Knee brake, Trendelenburg, Reverse Trendelenburg
Safe Working Load	573 lbs	Max Patient Weight	496 lbs


*Mattress platform extension infill piece MUST be used if the mattress platform is extended.



FloorBed® Accessories



Head and foot board bumpers



Mattress platform extension – increasing length to 84 inches



Assist bar to assist with bed mobility



Safety mat provides same height as a FloorBed® with 6 inch foam mattress



Mattress pump bracket

H&R Healthcare is passionate about service, striving to provide high quality products and superior customer care. Privately owned and operated for over 25 years, we have established ourselves as the premier provider of support surfaces, negative pressure wound therapy, and specialty medical equipment. From acute to long-term to home care settings, customers rely on H&R to provide cutting-edge healthcare products combined with clinical support, same-day delivery, and billing services.