


Gel Seat Cushion



The Gel Seat Cushion's unique design incorporates the strengths of gel, foam, and air, delivering outstanding pressure reduction, superior patient comfort, and proper patient positioning. The gel insert is sandwiched between multi-density foam layers to create a low interface pressure environment.

The positioning air bladder comes pre-filled from the factory and can be altered to specific patient/therapist preferences. The cover is manufactured to significantly reduce shear and friction.

The reticulated foam top layer is designed to aid in vapor transfer thus reducing tissue maceration. The T-Designed (3 cell) bladder allows proper gel displacement while providing optimum positioning. The cushion is a simple to use clinically superior pressure reducing/positioning seat cushion that will meet the needs of many patients.

Available Sizes: 18" x 18" , 24" x 17" , 30" x 19"

Additional sizes available upon request.