



Single Patient Disposable Fully Functioning Negative Pressure System

The UNO Pump

Lightweight - only 400 grams

Small - 7.6cm x 11.2cm x 6cm, easily and discreetly carried by the patient

Therapy Treatment Time

We offer two options with 7 or 15 days of effective therapy

2 Modes of Operation

1. **Continuous** - at either 125mmHg or 80mmHg
2. **Variable** - at either 80mmHg/30mmHg or 125mmHg/30mmHg

4 Alarm Settings

Leak, Canister Full, Blockage, Battery Low - you can be confident that the UNO system is providing optimal therapy at all times

Pump Carrying Case

- Washable & Reusable
- Adjustable waist or shoulder strap
- Provides additional discretion
- Wearable on the go

Indications

- Dehisced wounds
- Trauma Wounds
- Partial Thickness Burns
- Flaps and Grafts
- **Chronic Wounds** - Pressure Ulcers, Venous Ulcers, Diabetic Foot Ulcers



KIT Presentation

Each system comes with everything in one box! Additional dressings and canisters are available to order separately.

- **UNO 7** includes 1 Canister, 2 Dressing Kits, 1 Carry Pouch, 2 sets of batteries
- **UNO 15** includes 2 Canisters, 5 Dressing Kits, 1 Carry Pouch, 5 sets of batteries

The Port and Dressings

- **One luer-lock** connection to the pump
- **Foam Layer:** quick absorption and maintaining a moist environment at the wound site
- **Super Absorber Layer:** aids movement of exudate into the port
- **Silicone Wound Contact Layer:** soft & gentle application and removal, can be repositioned
- **Film Layer:** allows small amounts of air to pass through, transpires exudate away from the wound
- **Security Strips:** for use if dressing adherence is difficult
- **Pump Position:** works effectively positioned either above or below the wound site

Canister Presentation

70ml canister enables the system to easily manage exudate on both low and moderately exuding wounds. Each canister includes a gelling agent and deodorizing agent.

H&R Healthcare is passionate about service, striving to provide high quality products and superior customer care. Privately owned and operated for over 25 years, we have established ourselves as the premier provider of support surfaces, negative pressure wound therapy, and specialty medical equipment. From acute to long-term to home care settings, customers rely on H&R to provide cutting-edge healthcare products combined with clinical support, same-day delivery, and billing services.