

Powerful NPWT for treating wounds of all sizes for extended periods of time

The XLR8+ Pump

- **Lightweight** - only 585 grams, easy to maneuver and carry
- **Whisper-Quiet Operation** - patient gets better quality sleep
- **12-hour Lithium Battery** - easy to manage, minimal interruptions to care
- **Simple Lockable Controls** - no patient/caregiver identifiable interference with care

Advanced Variable Pressure

- Completely customizable high and low pressure settings
- Maintains negative pressure throughout the cycle, avoids losing airtight seal and displacement of the dressing
- Gentler on the patient and increases blood flow to the peri-wound
- In-growth of the granulation tissue to the foam dressing is virtually eliminated, making dressing changes less stressful and less painful to the patient

5 Alarm Settings

Leak, Canister Full, Blockage, Low Battery, Critical Battery

Green Foam

- **Reticulated Polyurethane Foam** - hydrolytic stability, increased performance
- **3-Dimensional Skeletal Strands** - enhanced filtering qualities
- **Non-Adherent Properties** - easier removal, almost eliminating pain
- Promotes perfusion and assists in tissue granulation
- Fraying is significantly reduced, preventing stray pieces from being left in wound bed

Silicone Port Pad

- **Softest Patient Interface** - enhanced patient comfort, less chance of pressure points
- **Film-Based Secure System** - easy application, dressing security
- **Integrated clips** - quick & easy changes, facilitates patient showering



Drape

- 3-Step Application Guide
- Pre-Marked Cutting Guides
- Clear Visibility
- Enhances Film Security




Canisters & Y-Connector

- **4 Capacity Canisters** - 400ml, 600ml, 800ml, 1100ml
- Slip and Clip attachment to pump
- Y-Connector allows multiple wound sites to be treated with one pump



Dressing Kits

- Thick Green Foam comes in various sizes - Small, Medium, Large, Extra-Large
- Spiro Foam Kit
- Bridge Dressing Kit



H&R Healthcare is a proud Medicare Part B Competitive Bid Winner

*Ability to go home with the patient, providing them with the **Full Continuum of Care***

H&R Healthcare is passionate about service, striving to provide high quality products and superior customer care. Privately owned and operated for over 25 years, we have established ourselves as the premier provider of support surfaces, negative pressure wound therapy, and specialty medical equipment. From acute to long-term to home care settings, customers rely on H&R to provide cutting-edge healthcare products combined with clinical support, same-day delivery, and billing services.